

## Product datasheet for **GA104816**

### Human TEAD3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | TEAD3 CRISPRa kit - CRISPR gene activation of human TEA domain transcription factor 3 |
| Symbol:       | TEAD3   |
| UniProt ID:   | <a href="#">Q99594</a>  |
| Synonyms:     | DTEF-1; ETFR-1; TEAD-3; TEAD5; TEF-5; TEF5  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)               |

#### Product images:



[View online >](#)