


Product datasheet for **GA104819**

Human TEF activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEF CRISPRa kit - CRISPR gene activation of human TEF transcription factor, PAR bZIP family member
Symbol:	TEF
Components:	GA104819G1, TEF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104819G2, TEF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104819G3, TEF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q10587
Synonyms:	KIAA1655
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)