

Product datasheet for **GA104821**

Human TIE2 (TEK) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEK CRISPRa kit - CRISPR gene activation of human TEK receptor tyrosine kinase
Symbol:	TEK
UniProt ID:	Q02763
Synonyms:	CD202B; TIE-2; TIE2; VMCM; VMCM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)