

Product datasheet for **GA104822**

Human TEP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEP1 CRISPRa kit - CRISPR gene activation of human telomerase associated protein 1
Symbol:	TEP1
Components:	GA104822G1, TEP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104822G2, TEP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104822G3, TEP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99973
Synonyms:	p240; TLP1; TP1; TROVE1; VAULT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)