

Product datasheet for **GA104825**

Human Telomerase reverse transcriptase (TERT) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TERT CRISPRa kit - CRISPR gene activation of human telomerase reverse transcriptase |
| Symbol: | TERT |
| UniProt ID: | O14746 |
| Synonyms: | CMM9; DKCA2; DKCB4; EST2; hEST2; hTRT; PFBMFT1; TCS1; TP2; TRT |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)