

Product datasheet for **GA104827**

Human Transferrin (TF) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TF CRISPRa kit - CRISPR gene activation of human transferrin
Symbol:	TF
UniProt ID:	P02787
Synonyms:	HEL-S-71p; PRO1557; PRO2086; TFQTL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)