

Product datasheet for **GA104833**

Human CP2c (TFCP2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TFCP2 CRISPRa kit - CRISPR gene activation of human transcription factor CP2
Symbol:	TFCP2
UniProt ID:	Q12800
Synonyms:	LBP1C; LSF; LSF1D; SEF; TFCP2C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)