

Product datasheet for **GA104843**

Human Transferrin Receptor 2 (TFR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TFR2 CRISPRa kit - CRISPR gene activation of human transferrin receptor 2
Symbol:	TFR2
UniProt ID:	Q9UP52
Synonyms:	HFE3; TFRC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)