

Product datasheet for **GA104844**

Human CD71 (TFRC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TFRC CRISPRa kit - CRISPR gene activation of human transferrin receptor
Symbol:	TFRC
UniProt ID:	P02786
Synonyms:	CD71; IMD46; p90; T9; TFR; TFR1; TR; TRFR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)