

Product datasheet for **GA104856**

Human TGF beta Receptor III (TGFB3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TGFB3 CRISPRa kit - CRISPR gene activation of human transforming growth factor beta receptor 3
Symbol:	TGFB3
UniProt ID:	Q03167
Synonyms:	betaglycan; BGCAN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)