

Product datasheet for **GA104860**

Human Transglutaminase 3 (TGM3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TGM3 CRISPRa kit - CRISPR gene activation of human transglutaminase 3
Symbol:	TGM3
UniProt ID:	<u>Q08188</u>
Synonyms:	TGE
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)