

Product datasheet for **GA104879**

Human TIMP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIMP1 CRISPRa kit - CRISPR gene activation of human TIMP metalloproteinase inhibitor 1
Symbol:	TIMP1
UniProt ID:	P01033
Synonyms:	CLGI; EPA; EPO; HCI; TIMP; TIMP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)