

Product datasheet for **GA104881**

Human TIMP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIMP3 CRISPRa kit - CRISPR gene activation of human TIMP metalloproteinase inhibitor 3
Symbol:	TIMP3
UniProt ID:	P35625
Synonyms:	HSMRK222; K222; K222TA2; SFD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)