

## Product datasheet for **GA104882**

### Human TIMP4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TIMP4 CRISPRa kit - CRISPR gene activation of human TIMP metalloproteinase inhibitor 4
Symbol:	TIMP4
UniProt ID:	<a href="#">Q99727</a>
Synonyms:	TIMP metalloproteinase inhibitor 4; tissue inhibitor of metalloproteinase 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)