

Product datasheet for **GA104889**

Human TLE 1 (TLE1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLE1 CRISPRa kit - CRISPR gene activation of human TLE family member 1, transcriptional corepressor
Symbol:	TLE1
UniProt ID:	Q04724
Synonyms:	ESG; ESG1; GRG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)