

Product datasheet for **GA104890**

Human TLE2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLE2 CRISPRa kit - CRISPR gene activation of human TLE family member 2, transcriptional corepressor
Symbol:	TLE2
UniProt ID:	Q04725
Synonyms:	ESG; ESG2; GRG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)