

Product datasheet for **GA104891**

Human TLE3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLE3 CRISPRa kit - CRISPR gene activation of human TLE family member 3, transcriptional corepressor
Symbol:	TLE3
UniProt ID:	Q04726
Synonyms:	ESG; ESG3; GRG3; HsT18976
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)