

Product datasheet for **GA104895**

Human Talin 1 (TLN1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLN1 CRISPRa kit - CRISPR gene activation of human talin 1
Symbol:	TLN1
UniProt ID:	Q9Y490
Synonyms:	ILWEQ; TLN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)