

## Product datasheet for **GA104899**

### Human TLR3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TLR3 CRISPRa kit - CRISPR gene activation of human toll like receptor 3
Symbol:	TLR3
UniProt ID:	<a href="#">O15455</a>
Synonyms:	CD283; IIAE2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)