

Product datasheet for **GA104900**

Human TLR4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLR4 CRISPRa kit - CRISPR gene activation of human toll like receptor 4
Symbol:	TLR4
UniProt ID:	O00206
Synonyms:	ARMD10; CD284; TLR-4; TOLL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)