

## Product datasheet for **GA104917**

### Human Tetranectin (CLEC3B) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CLEC3B CRISPRa kit - CRISPR gene activation of human C-type lectin domain family 3 member B
Symbol:	CLEC3B
UniProt ID:	<a href="#">P05452</a>
Synonyms:	TN; TNA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)