

Product datasheet for **GA104919**

Human Troponin C (TNNC2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNNC2 CRISPRa kit - CRISPR gene activation of human troponin C2, fast skeletal type
Symbol:	TNNC2
UniProt ID:	P02585
Synonyms:	fast skeletal muscle troponin C; OTTHUMP00000039531; troponin C2, fast; troponin C type 2 (fast)
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

