

Product datasheet for **GA104922**

Human TNFAIP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFAIP3 CRISPRa kit - CRISPR gene activation of human TNF alpha induced protein 3
Symbol:	TNFAIP3
UniProt ID:	P21580
Synonyms:	A20; AISBL; OTUD7C; TNFA1P2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)