

Product datasheet for GA104923

Human TSG6 (TNFAIP6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFAIP6 CRISPRa kit - CRISPR gene activation of human TNF alpha induced protein 6
Symbol:	TNFAIP6
UniProt ID:	P98066
Synonyms:	TSG-6; TSG6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)