

Product datasheet for **GA104924**

Human TNFRSF1A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF1A CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 1A
Symbol:	TNFRSF1A
UniProt ID:	P19438
Synonyms:	CD120a; FPF; MS5; p55; p55-R; p60; TBP1; TNF-R; TNF-R-I; TNF-R55; TNFAR; TNFR1; TNFR1-d2; TNFR55
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)