

Product datasheet for **GA104925**

Human TNFRSF1B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF1B CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 1B
Symbol:	TNFRSF1B
UniProt ID:	P20333
Synonyms:	CD120b; p75; p75TNFR; TBPII; TNF-R-II; TNF-R75; TNFBR; TNFR1B; TNFR2; TNFR80
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)