

## Product datasheet for **GA104926**

### Human Troponin C1 (TNNC1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNNC1 CRISPRa kit - CRISPR gene activation of human troponin C1, slow skeletal and cardiac type
Symbol:	TNNC1
UniProt ID:	<a href="#">P63316</a>
Synonyms:	CMD1Z; CMH13; TN-C; TNC; TNNC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)