

Product datasheet for **GA104927**

Human TNNI1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNNI1 CRISPRa kit - CRISPR gene activation of human troponin I1, slow skeletal type
Symbol:	TNNI1
UniProt ID:	P19237
Synonyms:	SSTNI; TNN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)