

Product datasheet for **GA104928**

Human Troponin I fast skeletal muscle (TNNI2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNNI2 CRISPRa kit - CRISPR gene activation of human troponin I2, fast skeletal type
Symbol:	TNNI2
UniProt ID:	<u>P48788</u>
Synonyms:	AMCD2B; DA2B; FSSV; fsTnI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)