

## Product datasheet for **GA104929**

### Human Cardiac Troponin I (TNNI3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNNI3 CRISPRa kit - CRISPR gene activation of human troponin I3, cardiac type
Symbol:	TNNI3
UniProt ID:	<a href="#">P19429</a>
Synonyms:	CMD1FF; CMD2A; CMH7; cTnI; RCM1; TNNC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)