

Product datasheet for **GA104930**

Human Troponin T1 (TNNT1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNNT1 CRISPRa kit - CRISPR gene activation of human troponin T1, slow skeletal type
Symbol:	TNNT1
UniProt ID:	P13805
Synonyms:	ANM; NEM5; STNT; TNT; TNTS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)