


Product datasheet for **GA104931**

Human Cardiac Troponin T (TNNT2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNNT2 CRISPRa kit - CRISPR gene activation of human troponin T2, cardiac type
Symbol:	TNNT2
UniProt ID:	P45379
Synonyms:	CMD1D; CMH2; CMPD2; cTnT; LVNC6; RCM3; TnTC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)