

Product datasheet for **GA104933**

Human Transition Protein 1 (TNP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNP1 CRISPRa kit - CRISPR gene activation of human transition protein 1
Symbol:	TNP1
Components:	GA104933G1, Transition Protein 1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104933G2, Transition Protein 1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104933G3, Transition Protein 1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P09430</u>
Synonyms:	TP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)