

## Product datasheet for **GA104935**

### Human TNR activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNR CRISPRa kit - CRISPR gene activation of human tenascin R
Symbol:	TNR
UniProt ID:	<u><a href="#">Q92752</a></u>
Synonyms:	TN-R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)