

Product datasheet for **GA104937**

Human TNXB activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TNXB CRISPRa kit - CRISPR gene activation of human tenascin XB |
| Symbol: | TNXB |
| UniProt ID: | P22105 |
| Synonyms: | EDS3; HXBL; TENX; TN-X; TNX; TNXB1; TNXB2; TNXBS; VUR8; XB; XBS |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)