

Product datasheet for **GA104941**

Human Topoisomerase III alpha (TOP3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TOP3A CRISPRa kit - CRISPR gene activation of human DNA topoisomerase III alpha
Symbol:	TOP3A
UniProt ID:	Q13472
Synonyms:	TOP3; ZGRF7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)