

Product datasheet for **GA104951**

Human Triosephosphate isomerase (TPI1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPI1 CRISPRa kit - CRISPR gene activation of human triosephosphate isomerase 1
Symbol:	TPI1
UniProt ID:	P60174
Synonyms:	HEL-S-49; TIM; TPI; TPID
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)