

## Product datasheet for **GA104954**

### Human Tropomyosin 3 (TPM3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TPM3 CRISPRa kit - CRISPR gene activation of human tropomyosin 3
Symbol:	TPM3
UniProt ID:	<a href="#">P06753</a>
Synonyms:	CAPM1; CFTD; HEL-189; HEL-S-82p; hscp30; NEM1; OK/SW-cl.5; TM-5; TM3; TM5; TM30; TM30nm; TPMsk3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)