

Product datasheet for **GA104959**

Human TPR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPR CRISPRa kit - CRISPR gene activation of human translocated promoter region, nuclear basket protein
Symbol:	TPR
UniProt ID:	<u>P12270</u>
Synonyms:	nuclear pore complex-associated protein TPR; nucleoprotein TPR; OTTHUMP00000033584; translocated promoter region (to activated MET oncogene); tumor potentiating region
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)