

Product datasheet for **GA104962**

Human TPTE activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPTE CRISPRa kit - CRISPR gene activation of human transmembrane phosphatase with tensin homology
Symbol:	TPTE
Components:	GA104962G1, TPTE gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104962G2, TPTE gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104962G3, TPTE gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56180</u>
Synonyms:	CT44; PTEN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)