

Product datasheet for **GA104969**

Human TRAF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRAF3 CRISPRa kit - CRISPR gene activation of human TNF receptor associated factor 3
Symbol:	TRAF3
UniProt ID:	Q13114
Synonyms:	CAP-1; CAP1; CD40bp; CRAF1; IIAE5; LAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)