

Product datasheet for **GA104979**

Human TRPC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPC3 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily C member 3
Symbol:	TRPC3
UniProt ID:	Q13507
Synonyms:	SCA41; TRP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)