

Product datasheet for **GA104980**

Human TRPC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPC4 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily C member 4
Symbol:	TRPC4
UniProt ID:	Q9UBN4
Synonyms:	HTRP-4; HTRP4; TRP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)