

Product datasheet for **GA104998**

Human TSSC3 (PHLDA2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHLDA2 CRISPRa kit - CRISPR gene activation of human pleckstrin homology like domain family A member 2
Symbol:	PHLDA2
UniProt ID:	<u>Q53GA4</u>
Synonyms:	BRW1C; BWR1C; HLDA2; IPL; TSSC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)