


Product datasheet for **GA105000**

Human GFUS activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSTA3 CRISPRa kit - CRISPR gene activation of human tissue specific transplantation antigen P35B
Symbol:	GFUS
UniProt ID:	Q13630
Synonyms:	FX; P35B; SDR4E1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)