

Product datasheet for **GA105003**

Human **TTC3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC3 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 3
Symbol:	TTC3
UniProt ID:	P53804
Synonyms:	DCRR1; RNF105; TPRDIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)