


Product datasheet for **GA105004**

Human **TTC4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC4 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 4
Symbol:	TTC4
Components:	GA105004G1, TTC4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105004G2, TTC4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105004G3, TTC4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O95801</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)