

Product datasheet for **GA105009**

Human TUB activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TUB CRISPRa kit - CRISPR gene activation of human TUB bipartite transcription factor
Symbol:	TUB
Components:	GA105009G1, TUB gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA105009G2, TUB gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA105009G3, TUB gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P50607</u>
Synonyms:	rd5; RDOB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)