

Product datasheet for **GA105022**

Human CD252 (TNFSF4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFSF4 CRISPRa kit - CRISPR gene activation of human TNF superfamily member 4
Symbol:	TNFSF4
UniProt ID:	P23510
Synonyms:	CD134L; CD252; GP34; OX-40L; OX40L; TXGP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)