

Product datasheet for **GA105023**

Human TNFRSF4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF4 CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 4
Symbol:	TNFRSF4
UniProt ID:	P43489
Synonyms:	ACT35; CD134; IMD16; OX40; TXGP1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)